

PBDE DETAILED FLOWCHART


1.




Remove upper rack from magnetic base. Label test tubes for Standards, Control, and Samples.

Tube #	Content
1, 2	Diluent/Zero Standard 0 ppb
3, 4	Standard 1 0.025 ppb
5, 6	Standard 2 0.05 ppb
7, 8	Standard 3 0.10 ppb
9,10	Standard 4 0.50 ppb
11,12	Standard 5 1.0 ppb
13	Control
14	Sample 1
15	Sample 2

Add 250 μ L of either Standards, Control or Samples to the bottom of each test tube by inserting the pipet tip all the way into the tube without touching the sides or the bottom of the tube.



7.




Do not separate upper rack from lower base. Using a smooth motion, *invert* the combined rack assembly over a sink and pour out the tube contents; keep inverted and gently blot the test tube rims on several layers of paper toweling.

8.



Add 1 mL of Washing Solution down the inside wall of each tube by using the technique described in Box 2. Wait 2 minutes. Using a smooth motion, invert the combined rack assembly over a sink and pour out the tube contents; keep inverted and gently blot the test tube rims on several layers of paper toweling. Repeat this step.



2.



Add 500 μ L of thoroughly mixed PBDE Antibody Coupled Magnetic Particles down the inside wall of each tube by using the technique described in Box 2. Vortex for 1 to 2 seconds (at low speed to minimize foaming).

3.



Incubate for 20 minutes at room temperature (15°-30°C).

4.




Add 250 μ L of PBDE Enzyme Conjugate down the inside wall of each tube by aiming the pipet tip 1/4" to 1/2" below the tube rim without touching the rim or tube wall with the pipet tip; deliver liquid gently.

9.



Lift the upper rack (with its tubes) off the magnetic base; add 500 μ L of Color Reagent down the inside wall of each tube by using the technique described in Box 2. Vortex for 1 to 2 seconds (at low speed to minimize foaming).



10.




React for 20 minutes at room temperature (15°-30° C). During this period, add 1 mL of Washing Solution into a clean tube for use as an instrument blank in Step 10.

5.



Incubate for 20 minutes at room temperature (15°-30°C).

6.




Combine the upper rack with the magnetic base; press all tubes into base; allow 2 minutes for the particles to separate.

11.



Add 500 μ L of Stopping Solution down the inside wall of each tube by using the technique previously described. Read results at 450 nm within 15 minutes after adding the Stopping Solution. Multiply results of extracted soil samples by the appropriate factor.









[Safety Caution: Stopping Solution contains 0.5% sulfuric acid.]

For Ordering or Technical Assistance Contact:
 ABRAXIS 54 Steamwhistle Drive, Ivyland, PA 18974
 Phone: 215-357-3911 Fax: 215-357-5232
 Website: www.abraxiskits.com E-mail: info@abraxiskits.com



PBDE Assay Kit Part # 500090 100 Tubes

PBDE CONCISE FLOWCHART

<p>1.</p>  <p>Separate the rack.</p> <p>Add 250 μL of either Standards, Control or Samples to the bottom of each test tube.</p>	<p>7.</p>  <p>Invert the combined rack.</p> <p>Blot.</p>
<p>2.</p>  <p>Mix the PBDE Antibody Coupled Paramagnetic Particles thoroughly and add 500 μl to each tube.</p> <p>Vortex for 1 to 2 seconds minimizing foaming.</p>	<p>8.</p>  <p>Add 1 mL of Washing Solution.</p> <p>Wait 2 minutes.</p> <p>Invert the combined rack.</p> <p>Blot.</p> <p>Repeat this step.</p>
<p>3.</p>  <p>Incubate for 20 minutes at room temperature.</p>	<p>9.</p>  <p>Separate the rack.</p> <p>Add 500 μL of Color Reagent to each test tube.</p> <p>Vortex.</p>
<p>4.</p>  <p>Add 250 μL of PBDE Enzyme Conjugate to each test tube.</p> <p>Vortex.</p>	<p>10.</p>  <p>Incubate for 20 minutes.</p> <p>Prepare blank.</p>
<p>5.</p>  <p>Incubate for 20 minutes at room temperature.</p>	<p>11.</p>  <p>Add 500 μL of Stopping Solution to each test tube.</p> <p>Read OD 450</p>
<p>6.</p>  <p>Combine the rack and magnetic base.</p> <p>Seat all tubes.</p> <p>Wait 2 minutes.</p>	

For Ordering or Technical Assistance Contact:
 ABRAXIS 54 Steamwhistle Drive, Ivyland, PA 18974
 Phone: 215-357-3911 Fax: 215-357-5232
 Website: www.abraxiskits.com E-mail: info@abraxiskits.com



PBDE Assay Kit Part # 500090 100 Tubes