

Algal Toxin Strip Test (Microcystins) Finished Drinking Water

1. Collect Sample




Collect 1 to 2 ml of sample.


2. Transfer

Using the pipette provided, add 7 drops of SAMPLE to the conical, flip-top tube.

(The conical, flip-top tube contains dried reagents.)



3. Shake and incubate



Close the conical, flip-top tube and shake for 30 seconds.



Incubate for 20 minutes.



(Dried reagents will dissolve, turning the sample purple.)

4. Test

Insert test strip into conical, flip-top tube with arrow pointing down.

Incubate for 10 minutes.




5. Dry

Remove test strip. Lay flat and allow to continue developing for 5 minutes.



6. Interpret



INTERPRET TEST

CONTROL LINE	TEST LINE	INTERPRETATION
NO CONTROL LINE PRESENT	NO TEST LINE PRESENT	INVALID RESULT
CONTROL LINE PRESENT	NO TEST LINE PRESENT	>5 ppb
CONTROL LINE PRESENT	MODERATE INTENSITY TEST LINE PRESENT	BETWEEN 0 AND 5 ppb

For Ordering or Technical Assistance Contact:
 ABRAXIS, LLC 54 Steamwhistle Drive, Warminster, PA 18974
 Phone: 215-357-3911 Fax: 215-357-5232 www.abraxiskits.com

5 Test: PN 520016
 20 Test: PN520017