

Importance of β -agonist Determination

Veterinary residues in foods pose a serious threat to public health. β -agonists are a group of veterinary drugs that have been used illegally in some countries. β -agonists act by impeding the uptake of adrenal hormones by nerve cells and through the stimulation of the cardiac system. They alter body composition by redistributing fat from muscle tissue, resulting in higher production efficiencies. The monitoring of raw meat and animal feed for drug and chemical residues is necessary to ascertain that these compounds are not misused and do not present a danger to consumers.


The potential risk for human health posed by the presence of β -agonists is high, due to the severity of the possible adverse effects. The β -agonist clenbuterol has been implicated in many poisoning cases in European and Asian countries. Although meat is the most frequently analyzed sample matrix, other matrices, such as water, milk and feed are also routinely analyzed.

The Abraxis β -agonist ELISA allows the determination of 41 water samples in duplicate determination. Only a few milliliters of sample are required. The test can be performed in approximately 1 hour.

Performance Data

Test sensitivity:

The limit of detection for Clenbuterol calculated as $X_n \pm 3SD$ ($n=16$) or as 90% B/Bind is equal to 0.065 ng/mL. LOQ is 0.10 ng/mL. The concentration of residue necessary to cause 50% inhibition (50% B/B₀) is approximately 1 ng/mL. Determinations closer to the middle of the calibration range of the test yield the most accurate results.



Test reproducibility: Coefficients of variation (CVs) for standards: <10%; CVs for samples: <15%.

Selectivity: This ELISA recognizes several β -agonists with varying degrees:

Cross-reactivities:	Clenbuterol	100%
	Bromobuterol	95%
	Cimaterol	95%
	Salbutamol	80%
	Ractopamine	80%
	Tulobuterol	50%
	Mabuterol	45%
	Zilpaterol	20%
	Phenethylamine A	115%

Samples: Water. See Preparation of Samples section.

General Limited Warranty: Abraxis, Inc. warrants the products manufactured by the Company against defects and workmanship when used in accordance with the applicable instructions for a period not to extend beyond the product's printed expiration date. **Abraxis makes no other warranty, expressed or implied. There is no warranty of merchantability or fitness for a particular purpose.**

For ordering or technical assistance contact:

Abraxis, Inc.
124 Railroad Drive
Warminster, PA 18974
Tel.: (215) 357-3911
Fax: (215) 357-5232
Email: info@abraxiskits.com
WEB: www.abraxiskits.com

R051117

β -agonist (BR) ELISA, Microtiter Plate

Enzyme-Linked Immunosorbent Assay for the Determination of β -agonists in Water Samples



Product No. 515535W

1. General Description

The β -agonist Broad Reactivity (BR) ELISA is an immunoassay for the detection of β -agonists in water samples. This test is suitable for the quantitative and/or qualitative detection of β -agonists in water samples (for other sample matrices use kits PN 515545 or PN 515535). Positive samples should be confirmed by HPLC, GC/MS, or other conventional methods.

2. Safety Instructions

The standard solutions in this test kit contain small amounts of Clenbuterol. In addition, the substrate solution contains tetramethylbenzidine and the stop solution contains diluted sulfuric acid. Avoid contact of stopping solution with skin and mucous membranes. If these reagents come in contact with the skin, wash with water.

3. Storage and Stability

The β -agonist ELISA should be stored in the refrigerator (4–8°C). The solutions must be allowed to reach room temperature (20–25°C) before use. Reagents may be used until the expiration date on the box.

4. Test Principle

The test is a direct competitive ELISA based on the recognition of β -agonist compounds by specific antibodies. β -agonists, when present in a sample, and a β -agonist-enzyme conjugate compete for the binding sites of anti- β -agonist antibodies immobilized on the wells of a microtiter plate. After a washing step and addition of the substrate solution, a color signal is produced. The intensity of the blue color is inversely proportional to the concentration of β -agonist present in the sample. The color reaction is stopped after a specified time and the color is evaluated using an ELISA reader. The concentrations of the samples are determined by interpolation using the standard curve constructed with each run.

5. Limitations of the β -agonist ELISA, Possible Test Interference

Numerous organic and inorganic compounds commonly found in samples have been tested and found not to interfere with this test. However, due to the high variability of compounds that might be found in samples, test interferences caused by matrix effects can not be completely excluded.

Mistakes in handling the test can also cause errors. Possible sources for such errors can be: Inadequate storage conditions of the test kit, incorrect pipetting sequence or inaccurate volumes of the reagents, too long or too short incubation times during the immune and/or substrate reaction, extreme temperatures during the test performance (lower than 10°C or higher than 30°C).

Each reagent is optimized for use in the Abraxis β -agonist ELISA kit (for water). Do not substitute reagents from any other manufacturer into the test kit. Do not combine reagents from other Abraxis β -agonist ELISA kits with different lot numbers.

The Abraxis β -agonist ELISA kit provides screening results. As with any analytical technique (GC/MS, HPLC, etc.), positive samples requiring some action should be confirmed by an alternative method.

Working Instructions

A. Materials Provided

1. Microtiter plate coated with an antibody against β -agonist.
2. Clenbuterol Standards (7): 0, 0.1, 0.25, 0.5, 1.0, 2.0, 8.0 ng/mL.
3. β -agonist-HRP Conjugate (250X concentrate), 0.1 mL/vial.
4. Conjugate Diluent Solution (ready to use), 16 mL
5. Wash Buffer (5X) Concentrate, 100 mL, must be diluted before use, see Test Preparation (Section C).
6. Color (Substrate) Solution (TMB), 12 mL.
7. Stop Solution, 12 mL.

B. Additional Materials (not included with the test kit)

1. Micro-pipettes with disposable plastic tips (10-200 and 200-1000 μ L)
2. Multi-channel pipette (10-250 μ L) or stepper pipette with plastic tips (10-250 μ L)
3. Container with 500 mL capacity (for diluted 1X Wash Buffer, see Test Preparation, Section C)
4. Deionized or distilled water
5. 4 mL glass vials with Teflon-lined screw cap
6. Microtiter plate washer (optional)
7. Microtiter plate reader (wave length 450 nm)
8. Shaker for microtiter plates (optional)
9. Centrifuge, capable of spinning at 2,500 x g or 1.2 μ m glass fiber filter (optional)
10. Timer
11. Tape or Parafilm

C. Test Preparation

Micro-pipetting equipment and pipette tips for pipetting the standards and the samples are necessary. We recommend using a multi-channel pipette or a stepping pipette for adding the conjugate, antibody, substrate and stop solutions in order to equalize the incubations periods of the solutions on the entire microtiter plate. Please use only the reagents and standards from one package lot in one test, as they have been adjusted in combination.

1. Adjust the microtiter plate and the reagents to room temperature before use.
2. Remove the number of microtiter plate strips required from the foil bag. The remaining strips are stored in the foil bag and zip-locked closed. Store the remaining kit in the refrigerator (4-8°C).
3. The standard solutions, enzyme conjugate, antibody, substrate and stop solutions are ready to use and do not require any further dilutions.
4. The β -agonist-HRP Conjugate provided concentrated. Before each assay, calculate the volume of conjugate needed Dilute only the amount necessary for the samples to be analyzed. Once diluted, the conjugate solution will only remain viable for one week (stored frozen). If additional samples are to be analyzed more than one week after reconstitution, a new vial of conjugate must be prepared. To reconstitute, add 24 μ L of HRP Conjugate to a vial containing 6.0 mL of Conjugate Diluent required and vortex thoroughly before use.
5. Dilute the Wash Buffer (5X) Concentrate at a ratio of 1:5. If using the entire bottle (100 mL), add to 400 mL of deionized or distilled water.
6. The stop solution should be handled with care as it contains diluted H₂SO₄.

D. Preparation of Water Samples

1. Transfer 1.0- 4.0 mL of water sample to a glass vial
2. Centrifuge glass vial for 5 minutes at 2500 x g, or filter samples using a 1.2 μ m glass fiber filter. If a filter or centrifuge is unavailable, allow the water sample to settle for at least 10 minutes prior to analysis.
3. The supernatant is ready to be analyzed in the assay (Assay Procedure, step 1).

The β -agonist concentration contained in the water samples is the same as the concentration obtained in ELISA (no dilution) Highly contaminated samples outside of the calibration range of the assay must be diluted further and re-analyzed. **Sensitivity in water = 0.1 ppb.**

For additional matrices please contact Abraxis, Inc. (Kit PN 515535 or 515545).

E. Working Scheme

The microtiter plate consists of 12 strips of 8 wells, which can be used individually for the test. The standards must be run with each test. Never use the values of standards which have been determined in a test performed previously.

	1	2	3	4	5	6	7	8	9	10	11	12
A	Std 0	Std 4	Sam 2									
B	Std 0	Std 4	Sam 2									
C	Std 1	Std 5	etc.									
D	Std 1	Std 5	etc.									
E	Std 2	Std 6										
F	Std 2	Std 6										
G	Std 3	Sam 1										
H	Std 3	Sam 1										

Std 0-Std 6: Standards

0; 0.1; 0.25; 0.50; 1.0; 2.0; 8.0 ppb

Sam1, Sam2, etc.: Samples

F. Assay Procedure

1. Add **30 μ L of the standards or water samples** into the wells of the test strips according to the working scheme given. We recommend using duplicates or triplicates.
2. Add **100 μ L of diluted enzyme conjugate** solution to the individual wells successively using a multi-channel pipette or a stepping pipette.
3. Incubate the strips for 30 minutes at room temperature.
4. After incubation, remove the covering and vigorously shake the contents of these wells into a sink. Wash the strips **three times** using the diluted wash buffer solution. Use at least a volume of **250 μ L of 1X wash buffer** for each well and each washing step. Remaining buffer in the wells should be removed by patting the inverted plate dry on a stack of paper towels.
5. Add **100 μ L of substrate (color) solution** to the wells. Cover the wells with parafilm or tape and mix the contents by moving the strip holder in a circular motion on the benchtop for about 30 seconds. Incubate the strips for 20 minutes at room temperature. Protect the strips from direct sunlight.
6. Add **100 μ L of stop solution** to the wells in the same sequence as for the substrate solution.
7. Read the absorbance at 450 nm using a microplate ELISA photometer within 15 minutes after the addition of the stopping solution.

G. Evaluation

The evaluation of the ELISA can be performed using commercial ELISA evaluation programs (4-Parameter (preferred) or Logit/Log). For manual evaluation, calculate the mean absorbance value for each of the standards. Calculate the %B/B₀ for each standard by dividing the mean absorbance value for each standard by the Zero Standard (Standard 0) mean absorbance. Construct a standard curve by plotting the %B/B₀ for each standard on the vertical linear (y) axis versus the corresponding β -agonist concentration on the horizontal logarithmic (x) axis on graph paper. %B/B₀ for samples will then yield levels in ppb of a β -agonist by interpolation using the standard curve. Samples showing lower concentrations of β -agonist compared to Standard 1 (0.25 ng/mL) are considered as negative. Samples showing a higher concentration than Standard 5 (4.0 ng/mL) must be diluted further to obtain accurate results.

Semi-quantitative results can be derived by simple comparison of the sample absorbances to the absorbances of the standards. Samples with lower absorbances than a standard will have concentrations of β -agonist greater than that standard. Samples which have higher absorbances than a standard will have concentrations of β -agonist less than that standard.

As with any analytical technique (GC, HPLC, etc.), positive results requiring regulatory action should be confirmed by an alternative method.