

# *Instruments & Accessories*

## *Colorimeter, 450 nm*



### **Features**

- ◇ Small, lightweight design. Weighs less than 200 grams. Suitable for field testing or on-the spot process monitoring.
- ◇ Battery operated. Four (4) AAA alkaline batteries; approximate life is 750 tests.
- ◇ The electronic and optical components are housed in a durable, impact-resistant case, and the optical system features a long-lasting LED light.
- ◇ Readings directly from standard 12 mm round reaction tubes (with supplied adapter), no need to purchase special cuvettes or deal with sample transfers.

**NOTE:** Colorimeter is manufactured by the Hach company.

### **Ordering**

**PN 475002      Colorimeter, 450 nm**

**Abraxis LLC**  
**54 Steamwhistle Drive**  
**Warminster, PA 18974**  
Phone: (215) 357-3911  
FAX: (215) 357-5232  
Email: [info@abraxiskits.com](mailto:info@abraxiskits.com)  
WEB: [www.abraxiskits.com](http://www.abraxiskits.com)

