

Instruments & Accessories

Microtiter Plate Strip Reader



Features

- ◇ **Automatic strip transportation.** 3 strip loading capacity or break-apart, single, and double strips. The carriers come in two styles, 3 X 12 for break-apart microwell strips, and an optional 3 X 8 for non-break apart microwell strips.
- ◇ **Bichromatic optics.** Standard model includes **four filters** (405, 450, 495, 620) or **six filters Visible model** (405, 450, 495, 620, 545, 600 nm) models and **Six filter UV model** (405, 450, 492, 630, 545, and 340 nm) are available for a nominal charge.
- ◇ Easy to use with step by step prompting and on board graphic printer.
- ◇ **On board software calculates single and multi-point calibration curves, regressions, and a variety of cut-off options.** Calculation modes are single point, point-to-point curve fit, linear regression (log, logit and linear), % absorbance, uptake and cut-off absorbance.
- ◇ Large non-volatile memory stores more than 50 user-entered tests.
- ◇ NRTL and CE certified. Available in both 115V and 230V AC. Switch located in the back to use at 110V (default setting) or 220V.
1-year Warranty

Ordering

PN 475007

Microtiter Plate Strip Reader

Abraxis LLC
54 Steamwhistle Drive
Warminster, PA 18974
Phone: (215) 357-3911
FAX: (215) 357-5232
Email: info@abraxiskits.com
WEB: www.abraxiskits.com

