



TUBE HOLDER

STEP 1




Pre-fold all scores so they are relaxed.

STEP 2



Snap 3 deck tabs into back panel so they 'click'.

STEP 3



Snap 2 side flaps over base tabs so they 'click'.

For ordering or technical assistance contact:

PN 472285

Abraxis LLC
54 Steamwhistle Drive
Warminster, PA 18974
Tel.: (215) 357-3911
Fax: (215) 357-5232
Email: info@abraxiskits.com
WEB: www.abraxiskits.com

